

BD SurePath™ Collection Vial

Important Training Information

New internal ledge enables easier removal of collection device head

- Insert head of broom-like device into the large opening in the collection vial
- Rotate the handle of the collection device and gently pull up to detach device head from handle
- **DO NOT** insert collection device head into the smaller opening in the vial



2-D barcode improves positive sample identification from collection to result

- After collecting your patient sample, remove the top barcode from the vial and attach to the requisition form
- The second barcode remains on the vial to link the sample to the patient information
- If your office is using electronic medical records, scan the barcode on the vial to automatically link patient data to the sample vial




Helping all people
live healthy lives

Easy, Effective Sample Collection


1. Collect

Collect the cytology sample using a broom-like device with detachable head. Follow manufacturer's collection instructions for detachable head device.



2. Drop

Use the interior rim of the new BD SurePath™ Collection vial to pull off the head of the broom-like device into the large opening in the collection vial.



3. Send

Place the cap on the vial and tighten. Remove the top barcode and place on the requisition form to ensure positive sample identification. Send the BD SurePath™ Collection Vial to the lab for processing.



BD Totalys™ System • MultiProcessor and SlidePrep

BD, BD Logo and all other trademarks are property of Becton, Dickinson and Company. ©2013 BD.

XEUR5235-13 Printed in Europe

BD Diagnostics
Tullastrasse 8-12
69126 Heidelberg
Germany
www.bd.com/europe/ds